

# ioGlutamatergic Neurons™

Human iPSC-derived  
glutamatergic neurons

Learn more about  
ioGlutamatergic Neurons

ioCells™



## About the cells

ioGlutamatergic Neurons have been precision reprogrammed from human induced pluripotent stem cells (iPSC) using opti-ox™ technology. Within days, cells convert consistently to mature, functional glutamatergic neurons characterised by >80% expression of glutamate transporter genes VGLUT1 and VGLUT2.


Glutamatergic neurons are delivered cryopreserved and ready-to-culture making them a high-quality human model for fundamental research, disease modelling and drug discovery.

## Benchtop benefits



### QUICK

Ready for experimentation as early as 2 days post-revival and form functional neuronal networks at 17 days.



### SCALABLE


Industrial scale quantities at a price point that allows the cells to be used from research to screening scale.



### EASY TO USE


Cells arrive programmed to rapidly mature upon revival. One medium is required in a two-phase protocol.

## Cells arrive ready to plate



ioGlutamatergic Neurons mature rapidly and form structural neuronal networks over 11 days (upper panel). 100X magnification.

Immunofluorescent staining on post-revival day 11 demonstrates homogenous expression of the pan-neuronal protein, MAP2 and glutamatergic neuron-specific transporter, VGLUT2 (lower panel).




ioGlutamatergic Neurons are highly characterised and defined, so you know exactly what is in every vial.

## Get reproducible results from every vial with high lot-to-lot consistency

### Whole transcriptome analysis demonstrates high lot-to-lot consistency across three manufactured lots

Bulk RNA-sequencing analysis was performed on three different lots of ioGlutamatergic Neurons on day 0, day 11 and day 18 post-revival. The experimental design included three operators, each handling one lot replicate of the different manufactured lots. (A) A principal component analysis (PCA) to assess gene expression variance between manufactured lots showed a tight clustering of the samples at each timepoint. (B) Differential expression test reveals no statistically significant differentially expressed genes across the three lots at day 11 ( $|\log_{2}FC| > 0.5$  and  $FDR < 0.01$ ).




**Differential Gene Expression at day 11**

Lot 2 vs Lot 3	0
Lot 1 vs Lot 2	0
Lot 1 vs Lot 3	0
iPSC vs Lot 1	5,735

## Robust and scalable cells suitable for high-throughput screening

### ioGlutamatergic Neurons show good suitability for high-throughput screening in 384-well format plates


Cytotoxicity CellTiter-Glo® (CTG) and TR-FRET (HTRF®) assays for AKT serine/threonine kinase 1 (AKT) and Huntingtin (HTT) proteins were performed on ioGlutamatergic Neurons in 384-well plates treated with tool compound (cmp) at day 9 post-revival. Compound titration results in a concentration response curve for all three assays (mean±sd of 2 replicates). The CTG assay shows an excellent average signal/background ratio and high suitability for HTS. The HTRF assays show a lower average signal/background ratio, due to lower assay sensitivity, indicating ioGlutamatergic Neurons are also suitable for HTRF assays. Data courtesy of Charles River Laboratories.



## ioGlutamatergic Neurons form the isogenic control for ioDisease Model Cells™. This isogenic pairing enables you to make true comparisons in your data, and confidently link genotype to phenotype.

### Comparison by microelectrode array (MEA) analysis of wild type ioGlutamatergic Neurons (WT) and Huntington's disease model ioGlutamatergic Neurons HTT 50CAG/WT (HD)

Map of the Firing Rate distribution for (A) WT and (B) HD; Network Firing Rate, recorded for 300 sec. for (C) WT, and (D) HD; 38 days in vitro. During development and maturation, cells in both cultures showed a gradual increase in spontaneous activity. Wild type neurons showed higher spontaneous activity than the disease model. The data demonstrate the value of pairing ioGlutamatergic Neurons with ioDisease Model Cells as an isogenic control to enable disease phenotype characterisation. Data courtesy of Charles River Laboratories.



## Product information

---

### Cat code

io1001

---

### Starting material

Human iPSC line

---

### Karyotype

Normal (46, XY)

---

### Seeding compatibility

6, 12, 24, 48, 96 & 384 well plates

---

### Shipping info

Dry ice

---

### Donor

Caucasian adult male  
(skin fibroblast)

---

### Vial size

Small: >1 x 10<sup>6</sup> viable cells  
Large: >5 x 10<sup>6</sup> viable cells

---

### Quality control

Sterility, protein expression (ICC)  
and gene expression (RT-qPCR)

---

### Differentiation method

opti-ox™ cell reprogramming

---

### Recommended seeding density

30,000 cells/cm<sup>2</sup>

---

### User storage

LN2 or -150°C

---

### Format

Cryopreserved cells

---

### Product use

ioCells™ are for  
research use only

---

### Applications

Drug discovery, neurotoxicology,  
high throughput screening,  
CRISPR Screening,  
3D bioprinting

## Who we are

bit.bio combines the concepts of cell programming and biology to provide human cells for research, drug discovery and cell therapy, enabling a new generation of medicines.

This is possible with our precision human cellular reprogramming technology opti-ox™ – a gene engineering approach that enables unlimited batches of any human cell to be manufactured consistently at scale

For general information,  
email [info@bit.bio](mailto:info@bit.bio)

To learn more,  
visit [www.bit.bio](http://www.bit.bio)