



## **CRISPR-Ready** **ioGlutamatergic Neurons™**

Human iPSC-derived  
glutamatergic neurons  
expressing Cas9 for rapid  
gene knockout generation

Learn more about  
CRISPR-Ready  
ioGlutamatergic Neurons



## About the cells

CRISPR-Ready ioGlutamatergic Neurons are precision reprogrammed iPSC-derived neurons containing a constitutively expressed Cas9 nuclease. This product has been designed for scientists looking to generate gene knockouts in a physiologically relevant human cell. These cells arrive ready for guide RNA (gRNA) delivery by day 1 post-thaw. Using our optimised lentivirus or lipid-based gRNA delivery protocol, users can maximise their knockout efficiency and start measuring readouts from gene knockouts and CRISPR screens within days.

## Benchtop benefits



**READY-TO-USE**

Defined and characterised human neurons constitutively expressing Cas9, ready for knockout experiments from day 1.



**QUICK AND EASY**

Generate readouts within days using a simple protocol for cell maturation and gRNA delivery.




**HIGH KNOCKOUT EFFICIENCY**

Optimised protocols for lipid or lentivirus based gRNA delivery ensure maximal knockout efficiency.

## Go from seeding to knockout to readout in days

CRISPR-Ready ioGlutamatergic Neurons arrive ready to use. With our optimised lentivirus or lipid-based guide delivery protocols, users can cut experimental timelines to a matter of days by no longer needing to invest time engineering and characterising their own Cas9 stable iPSC lines or optimising iPSC differentiation protocols.


With these cells, robust experimental readouts can be achieved by simply delivering gRNA against your target gene. Users do not require prior expertise in iPSC differentiation or gRNA delivery optimisation.



## High-efficiency single-gene knockouts can be achieved with both lentiviral and lipid-based gRNA delivery using CRISPR-ready ioGlutamatergic Neurons, observed by amplicon sequencing

### High-efficiency single-gene knockouts can be achieved with both lentiviral and lipid-based gRNA delivery using CRISPR-ready ioGlutamatergic Neurons, observed by amplicon sequencing


SOX11 indel formation was measured by amplicon sequencing in CRISPR-Ready ioGlutamatergic Neurons transfected or transduced with a gRNA targeting SOX11 delivered at either day 1 or day 3 after thawing. gRNAs targeting SOX11 were delivered by either lentiviral transduction or as a synthetic gRNA (10  $\mu$ M) using Invitrogen™ Lipofectamine™ RNAiMAX transfection reagent. DNA was harvested from the cells after 3 days of culture post transduction or transfection for amplicon sequencing of SOX11. Comparable knockout efficiencies were achieved for both lentiviral transduction and lipid-based transfection.




### High knockout efficiencies have been validated at the protein level by immunofluorescent staining

Immunofluorescent staining on CRISPR-Ready ioGlutamatergic Neurons transfected or transduced with a gRNA targeting SOX11 shows high percentage of knockout cells. gRNAs targeting SOX11 were delivered by either lentiviral transduction or as a synthetic gRNA (10  $\mu$ M) using Invitrogen Lipofectamine RNAiMAX transfection reagent at day 1 or day 3 (data not shown) after thawing. Cells were fixed for SOX11 immunofluorescence staining after three days of culture post-transduction or transfection. Comparable knockout efficiencies were achieved for both lentiviral transduction and lipid-based transfection. 100X magnification.

#### Lipid-based transfection, Day 1




#### Lentiviral transduction, Day 1



## Neurons at scale for CRISPR screens

### CRISPR-ready ioGlutamatergic neurons can be scaled to CRISPR screening workflows, and used to identify disease-relevant targets

A pooled single-cell CRISPR screen was designed to study the transcriptomic effects of knocking out 100 disease-relevant genes in CRISPR-ready ioGlutamatergic Neurons. Lentiviral gRNA transduction was carried out at day 3, and single cell gene expression analysis was performed at day 12. Cells that received non-targeting control gRNA (green dots) are evenly distributed across the UMAP, meaning no significant transcriptomic effects were observed. Cells that received gRNA targeting aminoacyl-tRNA synthetase (aaRSs) genes (purple dots) clustered together, indicating similar transcriptomic profiles. Pathway analysis of this cluster showed activation of the unfolded protein response, a pathway of therapeutic significance in Charcot-Marie-Tooth disease



## Product information

---

### Cat code

ioEA1090

---

### Starting material

Human iPSC line

---

### Karyotype

Normal (46, XY)

---

### Seeding compatibility

6, 12, 24, 48, 96 & 384 well plates

---

### Shipping info

Dry ice

---

### Donor

Caucasian adult male  
(skin fibroblast)

---

### Vial size

Small: >1 x 10<sup>6</sup> viable cells

---

### Quality control

Sterility, protein expression (ICC)  
and gene expression (RT-qPCR)

---

### Differentiation method

opti-ox™ cell reprogramming

---

### Recommended seeding density

30,000 cells/cm<sup>2</sup>

---

### User storage

LN2 or -150°C

---

### Format

Cryopreserved cells

---

### Product use

ioCells™ are for  
research use only

---

### Applications

Single gene knockouts,  
Combinatorial gene knockouts,  
Pooled CRISPR screens, Arrayed  
CRISPR screens, Drug discovery,  
High throughput screening

## Who we are

bit.bio combines the concepts of cell programming and biology to provide human cells for research, drug discovery and cell therapy, enabling a new generation of medicines.

This is possible with our precision human cellular reprogramming technology opti-ox™ – a gene engineering approach that enables unlimited batches of any human cell to be manufactured consistently at scale

For general information,  
email [info@bit.bio](mailto:info@bit.bio)

To learn more,  
visit [www.bit.bio](http://www.bit.bio)